



Ezi-Flow CSV Bottle Filling Operating Instructions

Follow these instructions for manually filling an
Ezi-Flow CSV Pharma Process Bottle



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2014



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2015



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2019

VERSION 1.0
DATE PRODUCED 04/04/22

1

- During this procedure
- you will need:

- 1 Pharma process bottle with flexible connector chute and Passive connection.
- 2 Ezi-Grip single use bag clamp
- 3 Red Tamper Evident Safety Tag
- 4 Filling Funnel – either stainless steel or polypropylene
- 5 CSV Passive Opening Tool
- 6 Pharma process bottle lifting clamp (Optional)

In addition, a number of mobile or benchtop filling stands are also available to make the Filling Process easier.





2:

First, locate the Pharma process bottle complete with Flexible Connector Chute and passive connection.



Note: Lifting aids are also available to help with the manual handling of the Pharma process bottle



3:

The Ezi-Flow CSV Passive has a Manual Interlock which prevents the Passive Slider from being removed. To remove the Passive Slider, simply depress the Passive Interlock Mechanism using the Passive Opening Tool. Now you can fully withdraw the Slider from the passive Body.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



Filling Stand

4:

Place Bottle into the filling stand. Passive square edge facing forward.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



5:

Fully insert the Filling Funnel through the Passive Body. The Filling Band must tightly wrap around the Filling Funnel and should be positioned towards the end of the Funnel Spout. This will create a “Clean Area” and prevent product from contaminating the Flexible Connector chute/bottle & Passive during the Filling Process.



6:

Product can now be carefully dispensed through the Filling Funnel into the bottle via the flexible connector chute.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



7:

When Filling is complete, carefully remove the Funnel, ensuring that no residual product is retained on the surface of the Funnel.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



Clean Area

Bag Clamp

8:

The Ezi-Grip Single Use Bag Clamp needs to be used to seal off the Flexible Connector Chute and prevent product from contaminating the Passive. Be sure to position the Bag Clamp in the “Clean Area” between the Passive and the Filling Band.

 Note: Ensure the clamp is orientated with the red toggle facing downwards and on the left-hand side. This is an important step and promotes best working practice for the discharge procedure.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



Red Safety Toggle number 4

9:

Carefully, push the two halves of the Bag Clamp together.

When they are fully engaged, slide the Red Safety Toggle, marked number 4, towards the centre of the clamp.


A close-up photograph of the EZI-DOCK Passive Slider mechanism. The slider is a white plastic component with a blue handle, positioned on a white track. A red interlock mechanism is visible on the left side of the track. The slider is currently in a retracted position. The background shows a white ceiling with a grid pattern and a blue wall.

Passive Slider

10:

The Passive Slider can now be re-inserted into the Passive Body. This will re-engage the Interlock.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



Red Tamper Evident Security Tag

11:

Finally, insert the Tamper Evident Security Tag into the holes in the Clamp Handle and pull tightly. This will prevent accidental release of the Bag Clamp.



12:

The filling procedure is now complete, and the Pharma Process bottle is now ready for discharge.

EZI-DOCK SYSTEMS LTD
Ezi-Flow CSV Bottle Filling



Ezi-Dock Systems Limited

Park Court, Park Lane Business Park, Kirkby in Ashfield,
Notts NG17 9GU United Kingdom

 Tel: **+ 44 (0) 1623 888 000**

 Email: **sales@ezidock.com**

 Web: **www.ezidock.com**