SISTO-C Komponenten

SISTO-C components

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Zubehör/Accessories

Leckanzeige
Leakage indicator

Endschalter
Limit switches

Verriegelung
Locking device

Zubehör/Accessories

Endschalter
Limit switches

Leckanzeige
Leakage indicator

Hubbegrenzung
Travel stop


Stellungsregler
Positioner

Magnetventil
Solenoid valve

Membranen / Diaphragms

EPDM

TFM/EPDM kaschiert/bonded

TFM/EPDM 2-teilig/2-piece

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SISTO-C

SISTO-CT

SISTO-CBAV

SISTO-CM

2.1 SISTO-C Membranen

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2.1 SISTO-C Diaphragms

SISTO's philosophy is to supply diaphragms which can be employed as universally as possible.

In practice, this means that we are able to cover a comprehensive spectrum of applications with just a few different diaphragm types. Different diaphragm types of the same size can be exchanged without replacing further components. Our proven diaphragms are of course subject to continuous development and improvement.

The unique SISTO design with enclosed diaphragm and spiral support ensures utmost operating reliability and durability.

The following diaphragm types are available for the SISTO-C valve range.

Elastomer diaphragms

SISTO soft elastomers are made from an EPDM rubber which has been designed to withstand an operating temperature of up to 140°C. The diaphragm material corresponds to FDA

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requirements, EU1935 and has been tested and certified in accordance with USP guidelines.

The diaphragms are well suited for use in “cold applications” such as in WFI systems, and thanks to their high ozone resistance they have proven extremely successful in CIP applications with all types of standard cleaning media. These EPDM diaphragms are also ideal for sterile steam applications.

A piece of fabric integrated into the EPDM rubber enables the diaphragm to be used under vacuum. The diaphragms can easily handle solids-loaden fluids and have a long service life.

Elastomer diaphragms bonded with TFM

Manufactured in a specially developed process, these composite diaphragms are made of an FDA compatible EPDM rubber which is firmly bonded to a TFM film featuring improved mechanical properties and lower permeability compared with conventional PTFE. This material quality also complies with FDA specifications and the EU1935.

The TFM's excellent chemical and thermal properties make this diaphragm type the first choice for applications with chemically aggressive media and high temperatures. This diaphragm type is also recommended for use in media which need to be protected against rubber abrasion and fluids containing oils and greases.

As the rubber is fully covered by the TFM sheet, even mechanical damage to the sheet does not result in cavity formation and therefore to undefined spaces which cannot be cleaned. This composite diaphragm is vacuum-resistant.

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2-piece TFM/EPDM diaphragms

Diaphragms with thicker TFM film are also available for SISTO-C valves. In order that the actuation force can be fully utilised without affecting the plastic diaphragm's long service life, the TFM diaphragm is backed up by a second diaphragm made of EPDM.



A metal pin embedded in the TFM diaphragm helps lifting the diaphragm and thus opens the valve. To avoid damage to the diaphragm during valve closure, the metal pin is fitted to the bonnet assembly by means of a floating nut.

Both TFM and EPDM comply with FDA specifications. This 2-piece TFM/EPDM diaphragm is perfectly suited to applications involving high chemical and thermal stress.

	Membranwerkstoffe Diaphragm materials	EPDM	TFM/EPDM kaschiert / bonded	TFM/EPDM 2-teilig / 2-piece
Konformität Conformity	USP	✓	✓	✓
	FDA	✓	✓	✓
	EU 1935	✓	✓	✓
Einsatzgrenzen Temperature range	max. Temperatur kurzfristig max. temperature short-term	140°C/284°F	140°C/284°F	160°C/320°F
	max. Temperatur langfristig max. temperature long-term	120°C/248°F	120°C/248°F	150°C/302°F
	Mindesttemperatur minimum temperature	-20°C/-4°F	-20°C/-4°F	-20°C/-4°F

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