



| REVISION HISTORY |             |      |          |
|------------------|-------------|------|----------|
| REV              | DESCRIPTION | DATE | APPROVED |
|                  |             |      |          |



SCALE  
1:2.5



**Replacement Diaphragm for :-**  
 OEM PART = 3.000-P2  
 USED ON 3" (DN80) Diaphragm Valve

|         | NAME       | DATE     |
|---------|------------|----------|
| DRAWN   | M. MC KEON | 08/21/14 |
| CHECKED |            |          |

**Material -**  
**EPDM BACKING**  
**TFM-DIAPHRAGM**

Complies with CFR 177.2600  
 and Conforms to  
 USP Class VI Biological Reactivity  
 Tests Section #87 and #88



|                           |                |              |
|---------------------------|----------------|--------------|
| TITLE<br><b>DD30-STFM</b> |                |              |
| SIZE<br>A                 | <b>ROMYNOX</b> | REV<br>C     |
| FILE NAME: DD30-STFM.dft  |                |              |
| SCALE: 1:2 (mm)           | ECO:0605       | SHEET 1 OF 1 |

Disclaimer - Diaphragm Direct® does not manufacture, represent, nor sell OEM valves or diaphragms manufactured by any other company. Diaphragm Direct® brand diaphragms are not sponsored by any other manufacturers. Any name, model, part number, brand, image, material or referral is for descriptive purposes only. Diaphragm Direct® brand replacement diaphragms are guaranteed to fit OEM equipment and meet or exceed OEM specifications. Diaphragm Direct® brand replacement diaphragms are covered by the standard Diaphragm Direct® limited product warranty.