LIFEMETATM MANUAL CRIMPING TOOL

Disposable high containment transfer system

CFD CFE

Lifemeta[™] manual crimping tools can realize aseptic disconnection of silicone or TPE tubing, quick disconnection of the tubing with the metallic sleeve into the edge die to prevent external microorganisms from entering the fluid path during disconnection.

Lifemeta[™] manual crimping tools provide aseptic disconnection in classified or non-classified environments. It protects the liquid inside the tubing from external conditions and seals the tubing without risk of leakage.



For full technical information and options
Or search CFD - CFE on www.romynox.nl



Features

- Fast disconnection time
- Seal and disconnect in less than 5 seconds
- Secure and trusted seal
- Testing validates that the aseptic seal is stronger than the tube
- Intuitive operation
- Disconnect in one step-just cut the metallic sleeve, reduces risk of error